## LOVIBOND MINIKIT

GB: Caustic Alkalinity P  
D: Hydroxid - Konzentration  
F: Concentration d'hydroxyde  
I: Concentrazione Idrossido  
E: Concentración de Hidróxido  
P: Concentração de Hidróxido  
NL: Hydroxide-concentratie  
DK: Hydroxid - konzentration

20 - 500 mg/l CaCO<sub>3</sub>

AF 415  
41 41 50

## GB

### Introduction

The Alkalinity P Minikit is for determination of caustic Alkalinity and it is based on the tablet count method. The test range is 20 - 500 mg/l CaCO<sub>3</sub>. The caustic (or hydroxide) Alkalinity is a measure of the water up to the pH-value 8.3 (colour change of phenolphthaleine), specifically due to hydroxides.

### Contents

calibrated sample container, 100 ml  
calibrated sample container, 30 ml  
Alkalinity P tablets, 100 pieces in foil  
Alkalinity P (BaCl<sub>2</sub>) tablets, 40 pieces in foil  
operating manual in 8 languages

### Instructions

Range 20 - 500 mg/l Sample Size 100 ml

1. Take a sample of the water under test in the sample container provided, filling to the 100 mark.
2. Add two Alkalinity P (BaCl<sub>2</sub>) tablets, stopper the container and shake until the tablets have disintegrated.  
The absence of a blue colour indicates zero caustic alkalinity and the test may be terminated at this stage.
3. If a blue colour is produced, add one Alkalinity P tablet to the solution, stopper the container and shake until the tablet has disintegrated. Continue adding Alkalinity P tablets one at a time in this manner until the colour of the solution changes from  
**BLUE to YELLOW.**
4. Note that the Alkalinity P (BaCl<sub>2</sub>) tablets are to condition the sample only and are not included in the number of tablets used for purpose of the calculation. Note the number of Alkalinity P tablets used. Multiply the number of tablets used by 40 and subtract 20 from the result. This gives the caustic alkalinity of the water expressed as mg/l CaCO<sub>3</sub>.

This calculation is represented by the formula below:

$$\text{Caustic Alkalinity as (CaCO}_3\text{)} = (\text{No. of Alkalinity P tablets} \times 40) - 20$$

### Cleaning

Thoroughly rinse out the sample container after use.

### Example

total tablets used: 2 tablets Alkalinity P (BaCl<sub>2</sub>) + 6 tablets Alkalinity P  
counted tablets: 6  
(6 x 40) - 20 = 220 mg/l CaCO<sub>3</sub>

### Note

The total alkalinity of water may be determined using the LOVIBOND Total Alkalinity Minikit AF 413. This test is applicable to the testing of natural waters, etc., but may also be used on boiler water samples if

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