Lovibond®

MINIKIT

- | | | | |
|-----------|---------------------------|-----------|-------------------------|
| GB | Total Hardness, 0 - 35 °e | D | Gesamthärte, 0 - 28 °dH |
| F | Dureté totale, 0 - 50 °f | I | Durezza totale |
| E | Durezza total | P | Durezza total |
| DK | Total hårdhed | NL | Totale hardheid |

0 – 500 mg/l CaCO₃

AF 424

41 42 40

GB Introduction

This Minikit is for determination of Total Hardness (as CaCO_3) in water. The test range is 0 - 500 mg/l CaCO_3 . Total Hardness in water is caused by divalent cations. Calcium and magnesium ions are usually the only once in significant concentrations.

Contents

calibrated sample container, 100 ml
calibrated sample container, 30 ml
Total Hardness tablets, 200 pieces in foil
operating manual in 8 languages

Instructions

Range	0 - 500 mg/l CaCO_3	Sample volume	50 ml
Range	0 - 250 mg/l CaCO_3	Sample volume	100 ml

1. Select the sample size 50 ml or 100 ml appropriate to the range under test. Rinse the container with the sample and take a sample of the correct size in it.
2. Add one Total Hardness tablet, stopper the container and shake the container until the tablet disintegrates. Continue adding Total Hardness tablets one at a time in this manner until the colour of the solution changes from **VIOLET** to **BLUE**.

3. Note the number of tablets used until the colour of the solution changed and calculate the result from the formula below appropriate to the sample volume taken:

Sample volume 50 ml: Total Hardness (mg/l CaCO_3) = (No. of Tablets x 40) - 20

Sample volume 100 ml: Total Hardness (mg/l CaCO_3) = (No. of Tablets x 20) - 10

Conversion Table

	Alkaline Earth Ions mmol/l	Alkaline Earth Ions m.equiv/l	mg/l CaCO_3	German Deg. $^{\circ}\text{dH}$	English Deg. $^{\circ}\text{e}$	French Deg. $^{\circ}\text{f}$
1 mg/l CaCO_3 =	0.01	0.02	1.00	0.056	0.07	0.10

Cleaning

Thoroughly rinse out the sample container and the stopper after each use.

Designation according to EC guidelines

Total Hardness – Order Code 51 51 61

The product is not subject to the labelling regulations according to the computation procedure of the "General EC Classification Directive for Preparations" in its most recent valid version.

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are Trademarks of the
The Tintometer Ltd.
Tintometer GmbH