# LOVIBOND

## MINIKIT

GB: Sulphate                    E: Sulfato  
D: Sulfat                        P: Sulfato  
F: Sulfate                       NL: Sulfaat  
I: Solfato                        DK: Sulfat

40 - 4000 mg/l

AF 431  
41 43 10

## GB

### Introduction

The test is based on the turbidity method using a single tablet reagent. The test range is 40 - 200 mg/l SO<sub>4</sub>. The range can be expanded to 4000 mg/l SO<sub>4</sub> if the sample will be diluted.

### Contents

double tube assembly  
stirring rod  
cleaning brush  
bottle with distilled water, 100 ml  
syringe, 5 ml  
Sulphate Turbidity tablets, 50 pieces in foil  
operating manual in 8 languages

### Assembled unit

Outer cylinder (sample)  
Inner cylinder (plunger)  
Fill level  
Turbidity scale relates to mg/l  
Black Dot Indicator

Take reading from shoulder

### Procedure

A: Range	Sample Size
40 - 200 mg/l	20 ml

1. Separate the double tube assembly.
2. Fill outer cylinder to the 20 ml mark with the water sample.
3. Add a Sulphate Turbidity tablet. Crush the tablet with a clean stirring rod and mix well.
4. Particles which settle at the bottom after the tablet has dissolved can be disregarded. A uniform cloudiness of the water sample indicates the presence of sulphate.
5. Insert the inner tube slowly into the outer cylinder and allow it to fill with the water sample through the two holes at the bottom.
6. Hold the outer cylinder up as far as possible (see illustration) to avoid light losses. The observer looks vertically from above the unit and moves the inner tube gently up and down until the black spot on the base just disappears. Hold the unit so that position of the inner tube is maintained.
7. Read the graduation mark on the outer tube corresponding to the position of the lower edge of the inner tube. If the edge is between graduation it is possible to make an estimate of the

concentration.

8. Multiply value with the factor 2. The result represents the sulphate level in mg/l SO<sub>4</sub>.

**B:**

Range	Sample Size	Factor multiply reading by
80 - 800 mg/l	5 ml	8
200 - 2000 mg/l	2 ml	20
400 - 4000 mg/l	1 ml	40

1. Separate the double tube assembly.
2. Take a sample of the water under test using the 5 ml syringe, filling to the appropriate volume selected (see above). Fill the content of the syringe into the outer cylinder.
3. Add distilled water from the bottle up to 20 ml.
4. Take the measurement as described before (point A 3 - 7).
5. Read the graduation mark and multiply this figure by the factor given in the table above. The result represents the sulphate level in the original sample expressed as mg/l sulphate SO<sub>4</sub>.

**Note**

After each measurement, thoroughly clean the outer and inner cylinder and stirring rod with the brush included in the Minikit.

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