# Lovibond® Water Testing Tintometer® Group

## SD 300 pH

### pH-Wert • pH Value

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Seite 3–27

**GB** Instruction Manual

Page 28–53

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declares that this product

Product name: **SD 300 pH**

that the device corresponds to the essential protection ratings established in the Regulations of the Council for the Approximation of Legislation for the member countries regarding electromagnetic compability (2004/108/EG) and the low voltage directives (2006/95/EG). The conformity to EMC are verified under observance of following standards.

**EN 61326-1: 2006** (table 3, class B)  
**EN 61326-1: 2006** (addendum A, class B)

This declaration is responsible for the manufacturer

Dortmund, January 20, 2013

---

Cay-Peter Voss, Managing Director

# Table of contents

1.	General information . . . . .	30
2.	Safety . . . . .	30
2.1	Intended use. . . . .	30
2.2	Safety signs and symbols . . . . .	30
2.3	Safety instructions . . . . .	31
3.	Product description . . . . .	32
3.1	Delivery Contents . . . . .	32
3.2	Operating and maintenance notices . . . . .	32
4.	Operating instructions . . . . .	32
4.1	Display elements. . . . .	32
4.2	Controls . . . . .	33
4.3	Connections. . . . .	34
4.4	Pop-up clip . . . . .	34
5.	Set-Up . . . . .	35
6.	Basics for measurement . . . . .	35
6.1	pH measurement . . . . .	35
6.2	Redox measurement. . . . .	36
6.3	rH measurement. . . . .	36
6.3.1	Manual pH value (and temperature) setting . . . . .	36
6.3.2	Automatic adoption of pH value from pH measurement (only with Logger deactivated) . . . . .	37
6.4	Calibration of the pH measurement . . . . .	37
6.4.1	The automatic temperature compensation for the calibration . . . . .	37
6.4.2	Performing the calibration . . . . .	38
7.	Configuration of the instrument . . . . .	40
8.	Data logger . . . . .	45
8.1	Manual recording ("Func-Stor") . . . . .	45
8.2	Automatic recording with adjustable "Func CYCL" cycle . . . . .	46
9.	Universal output . . . . .	47
9.1	Interface . . . . .	47
9.2	Analogue output . . . . .	48
10.	Adjustment of the instrument. . . . .	48
11.	GLP . . . . .	48
11.1	Calibration interval (C.Int). . . . .	48
11.2.	Calibration data memory (rEAd CAL) . . . . .	49
12.	Alarm („AL.“) . . . . .	49
13.	Real time clock („CLOC“) . . . . .	49
14.	Battery replacement . . . . .	50
15.	Error and system messages . . . . .	51
16.	Return and disposal . . . . .	52
16.1	Return . . . . .	52
16.2	Disposal . . . . .	52
17.	Technical data . . . . .	52

## 1. General information

Through this document carefully and familiarise yourself with the operation of the instrument before using it. Keep this document ready to hand and in the immediate vicinity of the instrument so that you or technical staff can refer to it at all times in case of doubt.

Assembly, set-up, operation, maintenance and shut-down may only be performed by technically qualified personnel. The technical personnel must carefully read and understand the operating manual prior to beginning all work.

The liability and warranty of the manufacturer for damages and consequential damages are voided in the event of improper use, non-observance of this operating manual, use by insufficiently qualified personnel as well as unauthorized changes to the instrument.

The manufacturer is not liable for costs or damages arising through the use of this instrument, especially in the case of improper use or misuse or faults to the connections or the instrument.

The manufacturer assumes no liability from printing errors.

## 2. Safety

### 2.1 Intended use

The instrument is designed for the measurement of pH and Redox Potentials using suitable electrodes. The electrode connection is made through a BNC socket.

Please note: Different electrode types are required for the pH and redox measurement. In addition, it is possible to connect a temperature sensor (Pt1000 or NTC 30k, with banana plug). The temperature measurement enables the automatic temperature compensation (ATC) of the pH, rH or mVH measurement and is displayed separately.

The safety instructions in this operating manual must be observed (see below).

The instrument may only be used under the conditions and for the purposes for which it was designed.

The instrument must be handled with care and used in accordance with the technical data (do not throw, strike, etc.). Protect the instrument from dirt.

### 2.2 Safety signs and symbols

Warning notices in this document are identified as follows:

	1. Warning! This symbol warns of imminent threatening danger where fatality, severe bodily injury and/or severe property damage may occur in the case of non-observance.
	2. Attention! This symbol warns of potential dangers or hazardous situations in which damage to the instrument and/or the environment may occur in the case of non-observance.
	3. Note! This symbol draws your attention to processes, which have a direct influence on the operation in case of non-observance or which can cause an unforeseen reaction.

## 2.3 Safety instructions

This instrument is built and tested in accordance with the safety provisions for electronic measurement instruments. The fault-free function and operational safety of the instrument can only be guaranteed if common, general safety precautions as well as the instrument-specific safety instructions in this operating manual are observed.

1. The function and operational safety of the instrument can only be adhered to under the climatic conditions specified in the chapter "Technical data". If the instrument is transported from a cold environment to a warm environment, a fault of the instrument function may arise due to the build-up of condensation. In this case, it is necessary to wait until the instrument's temperature adjusts to the room temperature before use.

2.  If it is suspected that the instrument can not be used without possibly imposing a danger, it should be turned off immediately and the potential danger be identified before the instrument is used again. The safety of the user may be diminished by the instrument if it

- exhibits visible damages.
- no longer works as specified.
- was stored for an extended period in unsuitable conditions.

In case of doubt, send the instrument to the manufacturer for repair or maintenance.

3. Be especially careful when connecting the circuitry to other devices the circuitry for connection to other instruments. Under certain circumstances, internal connections in foreign instruments (e.g. GND connection with earth) can lead to impermissible voltage potentials, which can impair the function of or even destroy the instrument itself or a connected instrument.

4.  This instrument is not suitable for safety applications, Emergency Stop equipment or applications in which a malfunction could cause injuries and material damage. If this notice is not observed, severe harm to the health and property damage may occur.

5.  This instrument may not be used in a potentially explosive environment. Operation in a potentially explosive environment causes an increased risk of detonation, fire or explosion as a result of spark formation.

### 3. Product description

#### 3.1 Delivery Contents

Included as standard:

- SD 300 pH sensor with 2 AAA-batteries
- Operating manual
- Additionally as a set in a case: Standard buffer solutions 90 ml each  
Measurement cell(s) pH/ temperature

#### 3.2 Operating and maintenance notices

1. Battery operation:

If "bAt" is shown in the lower display, the batteries are depleted and must be replaced. However, the instrument will still function for a short time. If "bAt" is shown in the upper display, the battery voltage is no longer sufficient for operation of the instrument and the battery is now fully depleted. For battery replacement, see the chapter "14. Battery replacement".

When storing the instrument in environmental temperatures above 50 °C, the battery must be removed. If the instrument is not used for an extended period, the battery should be removed. However, the time must be reset after when the instrument is turned back on.

2. The instrument and sensors/electrodes must be handled with care and used in accordance with the technical data (do not throw, strike, etc.). Plugs and sockets must be protected from dirt.

3. USB:

Make sure to only connect the USB interface cable to permissible components.


Operation with the USB 300 interface cable is recommended. If this configuration option is being used, the instrument is powered from the USB interface of the connected PC Or USB mains adapter.

4. Display values with cable break or no connected pH or redox electrode: If no electrode is connected or the connecting cable is defective, the mV or pH values are still shown. However, these are not valid measurement results!

### 4. Operating instructions

#### 4.1 Display elements

1. **Main display:** pH value,  
Redox value (mV, mVH),  
rH value

2. **Secondary display:** Temperature measurement value

3. Display arrows for **measurement value units**

4. Evaluation of the electrode or battery status

5. Display elements for representation of the minimum/maximum/saved measurement value


6. **atc arrow:** in operating mode shows "pH", "mVH" or "rH", depending on whether a temperature sensor is connected, and the **automatic temperature compensation** thereby activated.

7. **stab arrow:** indicates a stable measurement

8. **cal arrow:** indicates operating mode "pH", that means the instrument is in calibration process.

9. **logg arrow:** Logger is ready.  
Arrow blinks: Automatic recording (Logg CYCL) is active.

## 4.2 Controls

### On/Off switch, lamp

Press briefly: Activate lamp  
or switch on instrument

Press and hold: Switch off instrument

### set / menu:

Briefly press: in "pH", "rH" and "mVH" mode: manual temperature entry, if no temperature sensor is connect in "rH" mode, additional manual entry of the pH value  
press for 2 sec. (Menu): Call up of the configuration

### min / max:

Press briefly: Display of the minimum or maximum measured value

press for 2 sec.: Delete the respective value

### cal: only in "pH" mode:

Briefly press: Display of the electrode status (electrode symbol + bar display)

press for 2 sec.: Start the pH calibration


### store / enter:

Logger off: Hold and save the current measurement value ("HLD" in display)

Logger on: Operation of the logger – Data logger chapter

Set/Menu: Confirmation of entries, return to measurement

### 4.3 Connections

**Universal output:** Interface, supply, analogue output (see chapter)

**BNC socket:** Connection for pH or redox electrode. Waterproof to IP65 with appropriate cable!

**Banana socket:**  
Pt1000 connection or NTC 30k temperature sensor


**With electrodes with integrated temperature-sensor,** the banana socket is connected on the outside.

**With a separate reference electrode** it is connected on the inside.


### 4.4 Pop-up clip

**Handling:**


- Pull at label "open" in order to swing open the pop-up clip.
- Pull at label "open" again to swing open the pop-up clip further.

Pop-up clip closed


Pop-up clip at position 90°

Pop-up clip at position 180°

**Function:**


- The device with a closed pop-up clip can be plainly laid onto a table or attached to a belt, etc.
- The device with pop-up clip at position 90° can be set up on a table, etc.
- The device with pop-up clip at position 180° can be suspended from a screw or the magnetic holder GMH 1300.

Device attached to a belt

Device set up on a table


Device suspended from magnetic holder

## 5. Set-up

Connect the electrode, switch on the instrument with the  button.

After the segment test  the instrument briefly shows information about its configuration:

 If a zero point or slope adjustment was made  
(See chapter 7. Instrument Configuration.)

Remove the protective cap from the electrode. (Caution: The cap should contain KCL 3 M or a storage solution)

Then the instrument is ready for measurement.

## 6. Principles of Measurement

### 6.1 pH measurement

The pH value describes the acidic or alkaline behaviour of an aqueous solution.

pH values below 7 are acidic (the smaller the greater the acidity); values above 7 are alkaline; pH 7 = neutral.

From the negative decadic logarithm of the hydrogen ion activity (often this is approximately equal to the concentration of hydrogen ions) the following is calculated:

$$\text{pH value} = -\log_{10} \left( \frac{c(\text{H}^+) \cdot f(\text{H}^+)}{(1 \text{ mol/l})} \right) \quad \text{with} \quad \begin{array}{l} c(\text{H}^+): \text{Hydrogen ion concentration in mol/l} \\ f(\text{H}^+): \text{Activity coefficient of the hydrogen ions} \\ \text{(usually less than 1)} \end{array}$$

The abbreviation "pH" stands for pondus Hydrogenii (Latin pondus: "weight"; hydrogenium: "hydrogen").

In order to register the pH value of a solution, it should always be recorded together with the temperature measurement, e.g.:

pH 5.87; 22.8 °C.

Reason: The pH value of most liquids changes with the temperature.

The pH measurement is a very precise, but very sensitive measurement. The measured signals are very weak (high-ohmic), especially when measured in low ionic solutions. Therefore, it must be ensured that

- faults (electrostatic charging, etc.) are prevented,
- a stable measurement value is achieved through slow stirring,
- the plug contacts are kept clean and dry,
- electrodes (except for special, waterproof designs) are , as far as possible, not immersed below the shaft,
- the electrode is calibrated sufficiently often (see below) (The calibration frequency depends on the electrode and usage and may vary between every hour and several weeks.),
- an appropriate electrode is used.

## 6.2 Redox measurement

The Redox Potential (or: ORP) indicates the extent to which the measured sample has an oxidising or reducing effect with respect to the hydrogen standard electrode.

This potential is frequently utilised in swimming pools as a measurement for the disinfectant effect of a chlorination. The redox value is also an important parameter for aquaria, because fish can only survive within a specific redox range. The measurement also plays an important role in potable water preparation, wastewater monitoring and in industrial applications.

The measurement takes place with the diffused silver/silver chloride electrodes (reference system with 3 molar potassium chloride solution). It can be read directly (mV setting) or automatically with the mvH unit setting and converted with compensation for temperature to the "hydrogen standard electrode reference system".

Calibration comparable with the pH measurement does not take place with the redox measurement. However, the suitability of the electrodes can always be checked with redox test solutions.

## 6.3 rH measurement

The rH value is calculated from a pH and a redox measurement. It is used, for example, in order to describe the antioxidant effect of foods. This is a measurement for the capacity of foods to reduce harmful free radicals (bioelectronics according to Prof. Vincent).

In order to determine the rH value of your solution, proceed as follows:

### 6.3.1 Manual pH value (and temperature) setting

The measurements for pH and temperature (if no temperature sensor is connected) can be entered manually. To do so, briefly press the  button and enter the temperature value

with the  and  buttons. After pressing the button again the pH value can be

changed (see also manual setting of the temperature). The entry is confirmed with .

## 6.3.2 Automatic adoption of pH value from pH measurement (only with Logger deactivated)

Always make sure during the measurements that your pH and redox electrodes are in good condition and are clean and dry prior to their introduction to the solution.

First place the pH and redox electrodes and the temperature sensor in the solution and stir carefully.

### 1. Measuring the pH value:

Connect the pH electrode and the temperature sensor to the instrument.

Then set the instrument to pH measurement and carry out a calibration of the electrode as necessary (see 6.4 "Calibration of the pH measurement" and 7 "Configuration of the instrument").

Then measure the pH value of the solution and save the measurement with the "enter"

button. Do not switch off the instrument until the rH measurement is finished, otherwise the pH value will be deleted and must be entered manually.

### 2. Determining the rH value:

Now connect the redox electrode and configure the SD 300 pH for rH measurement.

Now the rH value of the solution appears in the main display and the secondary display alternates between the previously measured pH value and the temperature.

## 6.4 Calibration of the pH measurement

The electrode data of pH electrodes may be subject to large fluctuations due to ageing and individual variation. Therefore, the current calibration with buffer solutions must be checked prior to a measurement; in case of deviations, a new calibration must take place. (see also chapter 11 GLP)

Buffer solutions are liquids which have an exact pH value. For the calibration the following may be used:

- standard series ("set" Delivery Contents)

- DIN series (CAL dIn; pH 1.68 (A), pH 4.01 (C), pH 6.87 (D), pH 9.18 (F) and pH 12.45 (G))

- arbitrary buffer (CAL Edit; neutral buffer in the range 6.5 - 7.5pH).

The shelf life of buffer solutions is limited and is also significantly shortened due to insufficient rinsing and drying, among other things, when changing between solutions. This can result in erroneous calibrations! Therefore, use buffer solutions which are as fresh as possible for the calibration and rinse with deionised or distilled water!

### 6.4.1 The automatic temperature compensation for the calibration

Both the signal of the pH electrode and the pH buffer are temperature-dependent. If a temperature sensor is connected, the temperature influence on the electrode is automatically compensated for both the measurement and the calibration. Otherwise, the actual temperature of the respective buffer must be entered as precisely as possible (see below).

When working with the standard or DIN buffer series, the temperature influences of the buffers must also be compensated for. With a manual buffer selection, the pH values of the buffers are entered with the corresponding temperature in order to provide the most precise calibration possible.

## 6.4.2 Performing the calibration

**Please observe: A calibration can only be performed in the temperature range from 0 - 60 °C!**

If you have not done so yet, select the "pH" measurement function and activate the 1, 2 or 3-point calibration and the corresponding buffer series (Std, Edit or din), as necessary (see 7 "Configuration of the instrument").

Carefully pull off the protective cap from the electrode (Caution! Contains 3 M KCl!).

Rinse off the electrode with distilled water and dry.

**Start of the calibration: Press and hold the**

**button for two seconds.**

The command for the measurement of the 1st calibration solution appears in the display. The

calibration can be cancelled at any time with the

button. In this case, the previous calibration remains valid.

### 1. Calibration point 1: "Pt. 1"

Place the electrode and the temperature sensor (if present) in the neutral solution and carefully stir.

(Any arbitrary solution (e.g. pH 4) can be used for 1-point calibration)

As soon as a stable measurement has been determined, the instrument continues to the next point.

\*1)

**Without temperature sensor: Manual entry of Buffer 1 temperature**

Enter with the buffer temperature with buttons:

or

The value is adopted with

and the next calibration step is displayed.

The calibration is now finished with 1-point calibration; the electrode status is signalled with the left bar display.

### 2. Rinse the electrode with distilled or deionised water and dry

### 3. Calibration point 2: "Pt. 2" (only with 2 or 3-point calibration)

Place the electrode and the temperature sensor (if present) in the second buffer solution. (e.g. with standard series: pH 4.01 or pH 10.01) and carefully stir.

As soon as a stable measurement has been determined, the instrument continues to the next point.

\*1)

**Without temperature sensor: Manual entry of Buffer 2 temperature**

Enter with the buffer temperature with buttons:

or

The value is adopted with

and the next calibration step is displayed.

The calibration is now finished with 2-point calibration; the electrode status is signalled with the left bar display.

#### 4. Rinse the electrode in distilled or deionised water and dry

#### 5. Calibration point 3: "Pt. 3" (only with 3-point calibration)

Please observe that both an acidic and an alkaline calibration point are required for a 3-point calibration.

Place the electrode and the temperature sensor (if present) in the third buffer solution. (e.g. with standard series: pH 10.01) and carefully stir.

\*1) As soon as a stable measurement has been determined, the instrument continues to the next point.

**Without temperature sensor: Manual entry of Buffer 3 temperature**

Enter with the buffer temperature with buttons:

or

The value is adopted with  and the next calibration step is displayed.


The calibration is finished; the electrode status is signalled with the left bar display.

\*1) For manual buffer setting (CAL Edit), the pH value of the solution must be entered with

the  or  buttons. With solutions of the Standard and DIN series, the pH value of the respective solution is automatically recognised.

The value is adopted with  and the next calibration step is displayed.

## Error messages of the pH calibration:

	Neutral buffer is not permitted: - Electrode is defective - Incorrect buffer solution  - Defective buffer solution	Clean the electrode, recalibrate. If error occurs again -> Replace electrode Always use the neutral buffer as the first solution! (Exception: Single point calibration) Use a fresh buffer solution
	Transconductance is too low: - Defective electrode - Defective buffer solution	Replace electrode Use fresh buffer solutions
	Transconductance is too high: - Defective electrode - Defective buffer solution	Replace electrode Use fresh buffer solutions
	Incorrect calibration temperature	Calibration is only possible in the range from 0 - 60 °C

Permissible electrode data:

Asymmetry:  $\pm 55$  mV

Transconductance: -62 - -45 mV/pH

## 7. Configuration of the instrument

Some menu items are accessible depending on the current instrument setting (e.g. some are locked when the logger contains data).

To configure, press  and hold "menu"  for two seconds, to open the menu ("SEt" main display"). With "menu"  select the desired menu branch; with the  button you can toggle to the corresponding parameters, which you can then change (selection of parameters with ).

The adjustment of the parameters takes place with the  and  buttons.

Pressing "menu"  again switches back to the main menu and saves the settings. The configuration is completed with .



If the "menu" and "store" buttons are pressed together and held for longer than two seconds, the settings are reset to the factory settings. If data in the single value logger (logger: 'Func Stor') is first shown as 'rEAd Logg' menu: see also chapter 8 Data logger for this purpose. If no button is pressed for a period of more than two minutes, the configuration is cancelled. No changes made up to that point are saved!

Menu	Parameter	Values	Meaning		
		 			
	rEAd Logg: Read of the individual logger data, see Chapter 8!				
	<b>Set Configuration: General settings</b>				
			<b>Input: Selection of the measurement</b>		**
			"rH" arrow	Measurement of the rH value	
			"mV" arrow	Measurement of the mV value (REDOX and/or ORP)	
			"mVH" arrow	Measurement of the mV value based on the hydrogen system	
	"pH" arrow	Measurement of the pH value			
			<b>pH resolution: Resolution of the pH display</b>		
			0.1 - 0.001	One tenth pH - one thousandth pH	
			<b>Calibration: Selection of the number of calibration points</b>		
			1-Pt	1-point (only offset calibration, incline -59.2 mV/pH)	
			2-Pt	2-point (neutral + one additional buffer)	
			3-Pt	3-point (neutral + one acidic + one additional buffer)	
			<b>Calibration: Selection of the buffer series</b>		
			Std	Standard buffer series (pH 7, pH 4, pH 10) included in the "set" scope of delivery	
			din	DIN 19266 buffer series pH 1.68(A), pH 4.01(C), pH 6.87(D), pH 9.18(F), pH 12.45(G)	
			<b>Calibration: Selection of the buffer series</b>		
Edit			Arbitrary buffer, manual setting		
		<b>Calibration: Time interval for the calibration memory (factory setting: 30 days)</b>			
		1 - 365	Time interval for calibration memory (in days)		
		oFF	No calibration memory		
		<b>Unit t: Selection of the temperature unit</b>		**	
		°C:	All temperatures specified in degrees Celsius		
		°F:	All temperatures specified in degrees Fahrenheit		

Menu	Parameter	Values	Meaning		
	 				
	<b>Auto</b>	<b>Auto Hold: Automatic measurement processing (only effective with Logger = oFF mode)</b>			
		on	Automatic measurement processing (only with Logger = oFF) Auto Hold		
		oFF	Standard hold function at the push of a button (only with Logger = oFF)		
	<b>P.oFF</b>	<b>Auto Power-Off : Automatic instrument shut-off.</b>			
		1 - 120	Shut-off delay in minutes. If no button is pressed and no data traffic is taking place through the interface, the instrument shuts off automatically after the lapse of this time		
		oFF	Automatic shut-off deactivated (continuous operation)		
	<b>LiTE</b>	<b>Background lighting</b>			
		oFF:	No lighting		
		5 - 120	Automatic shut-off of lighting after 5 - 120 s (factory setting: 5 s)		
		on	Lighting always on		
	<b>Out</b>	<b>Universal output</b>			
		SEr	Serial interface activated		
		dAC	Analogue output activated		
		oFF	Interface and analogue output off -> minimal power consumption		
	<b>Adr.</b>	01.11 - 91	Base address of the unit for serial interface communication		
	<b>dA.C.0</b>	e.g. -2.00 - 14.00 pH	Input of the measurement at which the analogue output should output 0V, e.g. with 0.00 pH		
	<b>dA.C.1</b>	e.g. -2.00 - 14.00 pH	Input of the measurement at which the analogue output should output 1V, e.g. with 14.00 pH		

<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <b>SET</b>            Corr         </div>	<b>Set Corr: Adjustment of the measurements</b>		**		
	<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <sup>mV</sup>  <b>OFFS</b> </div>	<b>Zero correction/offset of the voltage measurement</b>		**	
		oFF	No zero correction of the voltage measurement		
		-10.0 - 10.0 mV	Zero correction of the voltage measurement in mV		
	<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <sup>mV</sup>  <b>SCAL</b> %         </div>	<b>Slope adjustment of voltage measurement</b>		**	
		oFF	No slope adjustment for voltage measurement		
- 5.000 - 5.000 %		Slope correction of voltage measurement in %			
<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <b>OFFS</b> °C         </div>	<b>Zero correction/offset of the temperature measurement</b>		**		
	oFF	No zero correction of the temperature measurement			
	-5.0 - 5.0%	Zero correction of the temperature measurement in °C			
<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <b>SCAL</b> °C %         </div>	<b>Slope adjustment of temperature measurement</b>		**		
	oFF	No slope adjustment for temperature measurement			
	-5.00 - 5.00%	Slope correction of temperature measurement in %			
<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <b>SET</b>            AL.         </div>	<b>Set Alarm: Adjustment of the alarm function</b>				
	<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <b>AL. 1</b> </div>	On / No.So	Measurement channel pH/mV/rH; Alarm on with sound / alarm on without sound		
		oFF	No alarm function for measurement channel pH/mV/rH		
	<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <b>A.LLo</b> </div>	e.g. - 2.00 - 14.00 pH	Min. alarm threshold pH/mV/rH (not with AL. 1. oFF)		
		e.g. -2.00 - 14.00 pH	Max. alarm threshold pH/mV/rH (not with AL. 1. oFF)		
	<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <b>AL. 2</b> </div>	On / No.So	Temperature measurement alarm on with sound / alarm on without sound		
		oFF	No alarm function for temperature measurement		
	<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <b>A.2Lo</b> </div>	- 10.0 - +110.0 °C	Min. alarm threshold temperature (not with AL. 2. oFF)		
<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <b>A.2Hi</b> </div>	- 10.0 - +110.0 °C	Max. alarm threshold temperature (not with AL. 2. oFF)			

	<b>Set Logger: Adjustment of the logger function</b>		**		
		<b>Selection of the logger function</b>		*	
		CYCL	<b>Cyclic:</b> Logger function of cyclical logger		
		Stor	<b>Store:</b> Logger function Single value logger		
		oFF	No logger function		
	0:01 - 60:00	Cycle time in [minutes:seconds] with cyclical logger	**		
	<b>Set Clock: Adjustment of the real time clock</b>				
		HH:MM	<b>Clock:</b> Adjustment of the clock hours:minutes		
		YYYY	<b>Year:</b> Adjustment of the calendar year		
		DD.MM	<b>Date:</b> Adjustment of the date Day.Month		
	<b>rEAd CAL: Read of the calibration data:</b> see Chapter 11.2 Calibration data memory (rEAd CAL)				

(\*) If there is data in the logger memory, parameters marked with (\*) cannot be opened. If they need to be changed, the data must first be deleted!

(\*\*) With the logger running, parameters marked with (\*\*) cannot be opened.

## 8. Data logger

The instrument has two different logger functions:

"Func-Stor": Manual measurement recording at the push of the "store" button.

In addition, a measurement entry (L-Id) is requested

"Func-CYCL": Automatic cyclical recording in the adjusted time interval  
The logger records two measurements per data set.

A data set is comprised of: pH, mV, mVH or rH measurement value  
Temperature measurement value  
L-ID measurement point (only with "Func-Stor")  
Time and date at which the data is saved

The software GSOF3050 (V3.0 or later) is required for the evaluation and transfer of the data; it can be used to start and adjust the logger function very easily.

With the logger function activated (Func Stor or Func CYCL), the hold function is not available; then the "store" button must be used for the logger operation.

### 8.1 Manual recording ("Func-Stor")

#### a) Record measurements manually:

When the "Func Stor" logger function is selected (see "Configuration of the instrument"), up to 1000 measurements can be saved manually:

 **Press briefly:** Data set is saved ("St. XX" is briefly displayed).  
XX is the number of the data set)

**Measurement point entry "L-Id":** Number of the measurement point

with  or  buttons. Select number from 0 - 19999.

The entry is confirmed by pressing store .

If the logger memory is full, the following appears 

#### b) Call up manual recording:

Saved data sets can be read in the GSOF3050 software and can also be observed in the instrument display itself.

 **Press and hold for two seconds:** The following appears in the display:

 **"rEAd LoGG" only appears if the data sets have already been saved! Without data sets, the configuration menu appears**

 **Press briefly:** Switch between measurement values, measurement point and date+time display of the data set

 or  Switch between the data sets

 Close the display of the data sets

### c) Delete manual recording:

If data has already been saved, it can be deleted with the store button:

**Press and hold for two seconds:** Open the delete menu

Switch the selection:

CLr  
no

Delete nothing (cancel the process)

CLr  
ALL

Delete all data sets

CLr  
LAST

Delete the most recently saved data set

Confirmation of the selection, end of the delete menu

## 8.2 Automatic recording with adjustable "Func CYCL" cycle

If the "Func CYCL" logger function was selected (see "Configuration of the instrument"), measurements are automatically recorded in the interval of the adjusted time cycle. The logger cycle time is adjustable from 1 s to 60 min (see "Configuration of the instrument"). Recordable data sets: 10000

### a) Start logger recording:

**Press and hold for two seconds:** Start selection, then again  recording is started automatically.

Each save process is signalled by a brief display of "St.XXXXX".

In this mode, XXXXX stands for the number of the data set. If the logger memory is full,

the recording is stopped automatically, the display shows **Lo66**  
**FULL**

### b) Stop logger recording:

**Press and hold for two seconds:** If a recording is in progress, the stop menu appears

Switch the selection:

StoP  
no

Do not stop the recording (cancel process)

StoP  
YES

Stop recording

Confirmation of the selection, end of the delete menu

If an attempt to turn off the instrument is made when a cyclical recording is running, a query of whether the recording should be stopped automatically appears. The instrument can only be switched off if the recording is stopped. The Auto Power Off function is deactivated when a recording is in progress!

### c) Delete logger recording:

**Press and hold for two seconds: If logger data is present and the recording was already stopped, the delete menu appears**

Switch the selection:

or

**CLr**  
no

Delete nothing (cancel the process)

**CLr**  
ALL

Delete all data sets

**CLr**  
LAST

Delete the most recently saved data set

Confirmation of the selection, end of the delete menu

## 9. Universal output

The output can be used either as a serial interface (for USB 300 interface adapter) or as an analogue output (0-1V). If the output is not required, it should be deactivated (Out OFF), because this function consumes the battery.

If the instrument is operated with the USB 300 universal interface adapter, the instrument is powered from this interface.

### Plug assignment:

4: external supply +5V, 50mA

3: GND

2: TxD/RxD (3.3V logic)

1: +UDAC, analogue output

Only suitable adapter cables are permissible (accessories)

## 9.1 Interface

With a USB 300 galvanically isolated interface converter (accessory), the instrument can be connected directly to a USB interface of a PC. The transmission takes place coded in binary format and is protected against transmission errors through extensive safety mechanisms (CRC). The following standard software package is available:

GSOFT3050: Operating and evaluation software for the integrated logger function

### The measurement instrument has two channels:

- Channel 1: Current value channel pH, mV or rH and base address
- Channel 2: Temperature value

The measurement/alarm/range values output through the interface are always shown in the adjusted display unit!

## 9.2 Analogue output

An analogue voltage of 0-1 V can be tapped at the universal output jack (Out dAC setting). With DAC.0 and DAC.1 the analogue output can be scaled very easily.

It must be ensured that the analogue output is not too heavily stressed, otherwise the output value can be falsified and the power consumption of the instrument increases sharply. Loads of up to approx. 10 kOhm are harmless.

If the display exceeds the value adjusted with DAC.1, 1V is output.

If the display undercuts the value adjusted with DAC.1, 0V is output.

In the event of an error (Err.1, Err.1, etc.), a voltage slightly above 1V is output at the analogue output.

## 10. Adjustment of the instrument

With offset and scale the measurement inputs can be adjusted, both the voltage measurement and the temperature measurement. Requirement: Reliable references are available (e.g. ice water, regulated precision water baths, etc.):

If an adjustment is made (deviation from factory setting), this is signalled with the message "Corr" when the instrument is switched on.

Default setting for offset and scale are 'off' = 0.0, i.e. inputs are not changed.

Zero point correction:

**Displayed value = measured value - offset**

Zero point and slope correction:

**Display =(measured value - OFFS) \* (1 + SCAL / 100)**

(Display °F = (measured value °F - 32°F - OFFS) \* (1 + SCAL / 100))

## 11. GLP

Regular monitoring of the instrument and accessories is a part of the GLP (Good Laboratory Practice). With pH measurements, the correct pH calibration, in particular, must be assured. To ensure this, the instrument assists you with the following described function.

Prerequisite for the use of the GLP functions is that the electrode is not replaced. The data is saved in the instrument, however, it is based on the respective electrode.

### 11.1 Calibration interval (C.Int)

You can specify a fixed interval at which the instrument automatically reminds that a new calibration should be performed and/or the calibration is no longer valid.

The length of the interval can be set dependant on your use and the stability of the electrode. As soon as the interval has lapsed, "CAL" blinks in the display.

## 11.2 Calibration data memory (rEAd CAL)


The last 16 calibrations with date and results are stored in the instrument and can be retrieved.

### Display the calibration data memory:

Saved calibration data can be read in the GSOF3050 software and can be observed in the instrument display itself:

**Press and hold for two seconds:**

The following appears in the display:  (configuration level)

**Press until the following appears:**  read cal. = "Read calibration data"

**Press briefly: Switch between**

- U.Asy = Asymmetry voltage in mV
- SL. 1 = slope acid in mV/pH \*1)
- SL. 2 = slope alkaline in mV/pH \*1)
- Date+time display of the data set

In parallel, the bar display shows the electrode evaluation of the corresponding calibration

or

Switch between the calibration data sets

Close the display of the calibration data sets

- \* 1) 1-point calibration: slope acid = slope alkaline = 59.16 mV/pH is assumed  
2-point calibration: slope acid = slope alkaline = determined slope  
3-point calibration: slope acid and slope alkaline are determined separately

## 12. Alarm ("AL.")

There are three possible settings:

Off (AL.oFF), On with sound (AL.on), On without sound (AL.no.So).

In the following cases an alarm is issued when the alarm function is active (on or no.So.):

- Lower alarm threshold (AL. Lo) undercut
- Upper alarm threshold (AL. Hi) exceeded.
- Sensor error
- Low battery (bAt)
- Err.7: System error (is always signalled with sound)

In the event of an alarm, the "PRIO" flag is marked in the instrument response with interface accesses.

## 13. Real time clock ("CLOC")

The real time clock is required for the chronological assignment of the logger data and the calibration time point. Therefore, check the settings as necessary.


## 14. Battery replacement

Before changing batteries, please read the following instruction and follow it step by step. Not following the instruction may cause harm to the instrument or the protection against ingress of water and dust may be lost!

Avoid unnecessary opening of the instrument!

1. Release the screws of the protective casing and remove it.
2. Open the 3 Phillips screws at the backside of the instrument.
3. Lay down the still closed instrument, so that the display side points upwards.  
The lower half of the housing incl. the electronics should be kept lying down during battery change. This avoids loss of the 3 sealing rings placed in the screw holes.
4. Lift off upper half of housing. Keep an eye on the six function keys, to be sure not to damage them.
5. Change carefully the two batteries (Type: AAA).
6. Check: Are the 3 sealing rings placed in lower half?  
Is the circumference seal of the upper half sound and clean?
7. Close the housing, taking care that it is positioned correctly, otherwise the sealing may be damaged. Afterwards press the two halves together, lay the instrument with display pointing downwards and screw it together again

**Take care to screw only until until you feel increasing resistance, higher screwing force does not result in higher water protection!**

## 15. Error and system messages

Display	Possible causes	Remedy
No display or confused symbols,  Instrument does not respond to keypad.	Battery is depleted	Insert new battery
	System error	Disconnect battery, wait a short time and put the batteries in again (Chapter 14)
	Instrument defective	Send in for repair
<b>Err.1</b>	Measurement range exceeded	Check: Is the measurement above the permissible measurement range of the sensor? -> Measurement value is too high!
	Sensor defect	Send in for repair
<b>Err.2</b>	Below Measurement Range	Check: Is the measurement below the permissible measurement range of the sensor? -> Measurement value is too low!
	Sensor defect	Send in for repair
<b>Err.7</b>	System error	Send in for repair
	Measurement range significantly exceeded or undercut	Measurement range significantly too high or too low
> CAL < CAL blinks in the upper display	The pre-adjusted calibration interval has lapsed or the last calibration was invalid	The instrument must be calibrated
	Neutral buffer not permissible	
	Incorrect buffer solution	Always use the neutral buffer as the first solution! (Exception: Single point calibration)
	Defective buffer solution	Use fresh buffer solution
	Defective electrode	Clean the electrode, re-calibrate. If error persists - > replace electrode
	Transconductance is too low	
	Defective electrode	Replace electrode
	Defective buffer solution	Use fresh buffer solutions
	Transconductance is too high	
	Defective electrode	Replace electrode
	Defective buffer solution	Use fresh buffer solutions
	Incorrect calibration temperature	Calibration is only possible in the range from 0 - 60 °C

If "bAt" blinks in the display, the battery is depleted. Measurement can continue for a brief time. If "bAt" remains in the display, the battery is fully depleted and must be replaced. No more measurements can take place.

## 16. Return and disposal

### 16.1 Return

All instruments which are sent back to the manufacturer must be free from measurement substance remains and/or other harmful substances. Measurement substance remains on the housing or on the sensor can endanger persons or the environment.

Use suitable transport packaging for the return of the instrument, especially if it is still a functioning instrument. Make sure that the instrument is protected with sufficient insulating material in the packaging.

### 16.2 Disposal

Dispose of the depleted batteries at a collection centre designated for this purpose. The instrument may not be disposed of with household waste. Should the instrument be disposed of, send it directly to us (with sufficient postage paid). We will dispose of the instrument properly and in an environmentally friendly manner.

## 17. Technical data

Measurement ranges	pH	-2.000 - 16.000 pH
	Redox / mV	-1999.9 - 1999.9 mV
		Based on hydrogen system: -1792 - +2207 mVH (at 25°C, DIN 38404)
	rH	0.0 - 70.0 rH
	Temperature	-10.0 - +110.0 °C, Pt1000 or NTC 30k 14.0 - 110.00 °F
Precision	pH	±0.005 pH
	Redox / mV	±0.05% FS
	Temperature	±0.2 K
Operating conditions		-25 to 50 °C; 0 to 95 % r.h. (non-condensing)
Storage temperature		-25 to 70 °C

Connectors	pH, redox	BNC socket, suitable for standard BNC and water-tight BNC cable add. connection for reference electrode: 4 mm banana jack
	Temperature	Pt1000 or NTC 30k over 4 mm banana jack Automatic sensor recognition
	Interface / ext. supply	4-pin socket for serial interface and supply, analogue output 0-1V
Input resistance	pH, redox	>1012 Ohm
Display		4 ½ digit, 7-segment status display for battery and electrode with bars, illuminated
pH calibration	Automatic	1, 2 or 3-point calibration, either DIN 19266 buffer or technical buffer
	Manual	1, 2 or 3-point calibration
GLP		16-calibration memory Adjustable calibration intervals (1 to 365 days, CAL warning after lapse)
Data logger		Real time clock Cyclical: 10,000 data sets, variable cycle: 1s - 60 min Single: 1000 data sets, with measurement point entry and date + time
Alarm		Sound/visual/interface
Add. functions		Min/Max/Hold
Housing		Robust PA6 GB30 housing, incl. protective sleeve
	Protection class	IP65, IP67
	Dimensions L*W*H [mm]	164 * 128 * 37 incl. protective sleeve, approx. 250 g incl. battery and protective sleeve
Power supply Power consumption		2*AAA battery, (included in Delivery Contents) 2.0 mA (with Out = Off, corresponding to 500 h), lighting ~10mA (switches off automatically)
Battery replacement display		Automatic with depleted battery "bAt"; "bAt" blinking = warning
Auto-Off function		If activated, the instrument switches off automatically when not used for an extended period of time (variable 1 - 120 min).
EMC		The instrument corresponds to the essential safety requirements, which are defined in the Council Directive for the harmonisation of the legal requirements of the Member States over electromagnetic compatibility (2004/108/EC). Additional error: <1%

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**Tintometer AG**
5212 Hausen AG